


Take Your Own Accurate Custom Impressions

- ✓ 510K & CE approved
- ✓ Completely copy teeth details
- ✓ High precision impression
- ✓ OEM & ODM available
- ✓ Easy to use




Good stability



Best flexibility and High tear resistance




Excellent Hydrophilic



Long term storage



Best recovery



High color contrast

PRODUCT PARAMETER

Recovery of elasticity:

>99.0%

Hardness:

58+/-2

Maximum affordable temperature:


200 degree centigrade

Max. compressive deformation:

<1.0%


APPLY TO

01




Clear aligner retainer

02




Clip on Snap on veneers

03



Custom mouth guard night guard

IMPRESSION MATERIAL USE FOR?



Dental impressions are used to create a close replica of either your set of teeth or your oral tissue or both. A plastic or metal "horseshoe" shaped tray which comfortably fits into your teeth and gums is chosen and inserted into your mouth. The tray can either be molded for the top, bottom or both teeth.

HOW TO USE?



Step 1



Step 2



Step 3



Step 4



Step 5



Step 6



Excellent dental impression




Shallow Dental Impression



The walls collapsed inward on the impression

PUTTY AVAILABLE WEIGHT OPTION


15g, 20g, 25g, 28g, 30g, 350g, 10Kg



OPTIONAL BASE COLOR




Catalyst Color: White



Base Color: Blue, Red, Orange, Green, Purple, Yellow OR Customization

DIFFERENT PACKAGES FOR PUTTY IMPRESSION MATERIALS



Jelly cup package



White jar



Clear jar



Semi-up

OEM SERVICE



A. Custom Logo On The Box

- UV, Stress Convex Printing
- Hot Stamping



B. Custom Logo On The Jar

- Double Color Printing
- Custom Sticker Available



C. User Manual

- Language Customization
- Shapes, Color



D. Impression Tray

- Logo Printing
- Different Color Tray
- Small, Middle, Large Size



E. Putty Customization

- Setting Time, Mixing Time
- Biting Time
- Hardness, Color OEM



F. Extra Our Service

- Free Design Service
- USA Clearance